

Stealth Health AI Agent – Project Submission

Candidate

Shaswat

1. Project Overview

Stealth Health AI Agent is an AI-powered wellness coaching assistant designed to personalize health guidance using structured patient onboarding, adaptive daily check-ins, session memory, and protocol-grounded responses.

The system accepts unstructured patient onboarding information, extracts structured wellness data, adapts conversations based on the user's protocol timeline, and answers wellness questions strictly using a predefined protocol document.

2. Live Application

Live Streamlit Application

<https://stealth-health-ai-agent-qv84dd6fsmh4wrhureqd4z.streamlit.app/>

Query Parameter Demo

The application supports URL-based onboarding.

Example:

https://your-app-name.streamlit.app/?raw_data=I am Guest , 22 years old. My goal is improving sleep. I currently sleep 5 hours. Put me on Day 5.

Expected Extraction Result

```
{  
  "name": "Guest",  
  "age": 22,
```

```
"primary_goal": "improving sleep",  
"sleep_hours": 5,  
"current_day": 5  
}
```

3. GitHub Repository

Repository URL:

<https://github.com/shaswatgithub/Stealth-Health-AI-Agent>

The repository contains:

- Source Code
 - Agent Logic
 - Streamlit Application
 - Deployment Configuration
 - Documentation
 - Requirements File
-

4. Problem Statement

Build a Health Coach AI Agent that:

- Accepts patient onboarding data
 - Extracts structured patient information
 - Performs adaptive daily check-ins
 - Maintains session memory
 - Answers protocol questions with minimal hallucination
 - Uses a warm, clear coaching style
-

5. Solution Architecture

Workflow

User Input

↓

Gemini Extraction Agent
↓
Structured Patient Profile
↓
Session Memory
↓
Day-Based Coaching Logic
↓
Protocol Grounding
↓
Personalized Response

6. Key Features

A. Structured Patient Onboarding

The system extracts:

- Name
- Age
- Wellness Goal
- Sleep Habits
- Protocol Day

from natural language onboarding text.

Example Input:

"I am Shaswat, 22 years old. My goal is improving sleep. I currently sleep 5 hours. Put me on Day 5."

B. Adaptive Daily Check-ins

The AI dynamically changes behavior based on protocol progression.

Day 1

- Welcome message
- Goal review
- Habit baseline setup

Day 2–4

- Progress tracking
- Habit consistency monitoring

Day 5+

- Accountability coaching
 - Goal-focused follow-up
 - Obstacle identification
-

C. Session Memory

The application stores previous messages using Streamlit Session State.

This enables:

- Context retention
 - Personalized conversations
 - Multi-turn coaching interactions
-

D. Protocol-Grounded Responses

The AI is restricted to responding only from the provided wellness protocol.

If information is unavailable in the protocol, the assistant responds:

"I couldn't find that in the protocol document."

This significantly reduces hallucination risk.

7. Technology Stack

Component	Technology
Frontend	Streamlit
AI Model	Gemini 2.5 Flash
Data Validation	Pydantic
Session Memory	Streamlit Session State

Deployment	Streamlit Community Cloud
Source Control	GitHub

8. Agent Design

Agent State

The system maintains a structured patient profile:

```
{  
  "name": "Guest",  
  "age": 22,  
  "primary_goal": "Improve Sleep",  
  "sleep_hours": 5,  
  "current_day": 5  
}
```

Memory

Previous interactions are stored in session state and supplied back to the model for contextual continuity.

Reasoning

The agent combines:

- Patient Profile
- Timeline Context
- Protocol Rules
- Conversation History

to generate personalized responses.

9. Technical Decisions

Why Streamlit?

- Fast development
- Built-in chat interface
- Easy deployment
- Excellent demo capability

Why Gemini 2.5 Flash?

- Fast inference
- Strong instruction following
- Cost-efficient
- Reliable structured extraction

Why Pydantic?

- Schema validation
- Type safety
- Clean structured system state

Why Session State?

- Lightweight memory
 - No external database required
 - Suitable for MVP requirements
-

10. Future Improvements

Potential future enhancements include:

- PDF-based RAG pipeline
 - Long-term memory storage
 - User authentication
 - Progress analytics dashboard
 - Multi-user support
 - Vector database integration
 - Habit tracking visualizations
-

11. Demo Video

Video Walkthrough:

[PASTE VIDEO LINK HERE]

Suggested Demo Flow:

1. Open application
 2. Demonstrate onboarding extraction
 3. Show parsed patient profile
 4. Demonstrate Day 1 behavior
 5. Demonstrate Day 5 behavior
 6. Show protocol-grounded responses
 7. Demonstrate session memory
-

12. Conclusion

Stealth Health AI Agent successfully satisfies the MVP requirements by combining structured information extraction, adaptive coaching behavior, conversational memory, and protocol-grounded responses within a simple, deployable Streamlit application.

The solution demonstrates a practical implementation of an AI agent architecture for personalized wellness coaching while maintaining transparency and minimizing hallucinations.